

Cracking the Code on Nitrogen Methods

Sam Johnson
Inorganic Chemist

Why Nitrogen Testing in Wastewater Is Challenging (and Critical)

- Nitrogen isn't "one test" – it's a family of measurements
 - Ammonia Nitrogen (NH_3 as N)
 - Nitrate + Nitrite Nitrogen (NO_3/NO_2 as N)
 - TKN (Total Kjeldahl Nitrogen)
 - Total Nitrogen (TN)
- Environmental impacts
 - Eutrophication and Hypoxia/Anoxia
- Human health risks
 - Nitrates interfere with the body's ability to transport oxygen
 - Excessive nitrogen fertilizer contributes to soil nutrient depletion

Products Highlighted (Core Nitrogen Focus)

Primary Nitrogen Parameters


- Simple Nutrients (PT # 584 / QC # 505)
 - NH_3 as N (1 – 20 mg/L)
 - NO_3 as N (2 – 25 mg/L)
 - $\text{NO}_3 + \text{NO}_2$ as N (2.5 – 25 mg/L)
 - ortho PO_4 as P (0.5 – 5 mg/L)
 - Total Nitrogen (3 – 45 mg/L)
- Complex Nutrients (PT # 579 / QC # 525)
 - TKN (3 – 35 mg/L)
 - Total Phosphorus as P (0.5 – 10 mg/L)
- Nitrite (PT # 888 / QC # 770)
 - NO_2 as N (0.4 – 4 mg/L)
- Available as QR if needed



Understanding Nitrogen Forms in Wastewater

Nitrogen Cycle

- Ammonia ($\text{NH}_3/\text{NH}_4^+$)
- Nitrate (NO_3^-)
- Nitrite (NO_2^-)
- Biological nitrification and denitrification




Overview of Common Nitrogen Test Methods

Ammonia / Ammonium

- SM 4500-NH3 (Ion Sensitive Electrode)
 - Ion Strength Adjustment (ISA) pH to > 11 by converting NH_4^+ to NH_3
 - No distillation needed
 - pH electrode inside sensor
 - Gas permeable membrane
 - Filling solution reacts with NH_3 to produce pH shift
- Low concentrations and drift of mV readings most common challenges

Waters™ | ERA



Ammonia Methods Comparison

Nessler Reagent

- Colorimetric; reaction of ammonia with Nessler's reagent ($K_2[HgI_4]$) in an alkaline solution to form a yellow-to brown complex.
- 0.05–1.0 mg/L NH_3 -N
- High sensitivity, standardized (EPA Method 350.2), relatively straightforward post distillation.
- Contains toxic mercury, potential for interferences from amines and other compounds, distillation can be time-consuming.

Berthelot (Indophenol)

- Colorimetric; reaction of ammonia with phenol and hypochlorite in an alkaline solution to form a blue-colored indophenol.
- 0.025–5.0 mg/L NH_3 -N
- High sensitivity and accuracy, less toxic reagents than Nessler's method.
- Longer reaction time (1-2 hours), sensitive to pH and interferants.

Ion-Sensitive Electrode

- Potentiometric; measurement of the potential difference across an ammonia specific membrane.
- Wide range, typically 0.01 to 17,000 mg/L NH_3 -N.
- Wide dynamic range, real-time measurements, portable
- Potential for interference from volatile amines, requires careful frequent calibration and temperature control to maintain precision.

Ammonia Methods Comparison

- Nessler USEPA Method 8038 (Hach kits)
 - UV Vis spectrophotometry
- Berthelot (Indophenol)
- ISE
- Ion Chromatography
 - Cations column
 - Automated
 - No hazardous waste



Waters™ | ERA



Nitrate and Nitrite


- Nitrate is a regulated contaminant with a Maximum Contaminant Level (MCL) of 10 mg/L in drinking water
- EPA 300.0
- EPA 9056A
- EPA 353.2 (Colorimetry)
- SM 4500 NO2-B
- SM 4500 NO3-D (ISE)
- SM 4500 NO3-E (Manual Cd reduction)
- SM 4500 NO3-F (Automated Cd reduction)
- Hach 10206 and 10207



Total Kjeldahl Nitrogen (TKN)

- Represents the sum of organic and ammonia nitrogen in a sample, excluding nitrate and nitrite
- EPA 351.2 / SM 4500-N
 - Sulfuric Acid / copper and potassium sulfate
 - Heat to 380°C to convert to ammonium sulfate
 - Run on ISE
- Hach Method 10242
 - Digest samples at 120°C
 - Colorimetric
 - Also allows for total nitrogen

Waters™ | ERA



Data Reporting

Required reporting “as N”

- Report in mg/L units
- Use “as N” calibration standards
- H – 1.0079 g/mol
- N – 14.0067 g/mol
- O – 15.9994 g/mol
- NH_3 – 17.0304 g/mol ($\text{N} = \text{NH}_3 / 1.2159$)
- NH_4 – 18.0383 g/mol ($\text{N} = \text{NH}_4 / 1.2878$)
- NO_2 – 46.0055 g/mol ($\text{N} = \text{NO}_2 / 3.2845$)
- NO_3 – 62.0049 g/mol ($\text{N} = \text{NO}_3 / 4.4268$)
- Use CRMs to confirm your process is reading correctly

Waters™ | ERA



Emerging Trends and Best Practices

- When submitting PT's, use the name and analyte code stated on your accreditation certificate or permit
- Know what the concentration ranges are for your PTs and if the method can achieve them.
- Clear understanding of nitrogen species in wastewater
- Confidence in selecting appropriate test methods
- Improved data quality and regulatory confidence
- Practical guidance for day-to-day lab and plant decisions
- Monitor your calibration and QC data – great indicators of potential quality issues
- Run CRMs alongside PTs

An abstract graphic featuring a network of interconnected nodes and lines, resembling a molecular structure or a data network, set against a blue gradient background. The nodes are represented by small circles of varying sizes, and the lines are thin, light blue. The overall composition is centered and extends across the width of the image.

Waters™ |  ERA™

Thank you!