


Getting pH Right


Christy Abbas, PhD
Sr Product Marketing Manager

Outline


- What is pH and why do we care
- pH proficiency testing
- Waters ERA pH products
- pH analysis
 - Measurement
 - Calibration
 - Probe care and troubleshooting



Hydronium Ion (H⁺)




Hydroxide ion (OH⁻)



What is pH?

- pH denotes the “potential of hydrogen”
 - **p** is mathematical notation for $-\log$
 - **H** is the elemental symbol for Hydrogen
 - $\text{pH} = -\log[\text{H}^+]$
 - Concentration of H^+ is expressed in molar
- $10,000 = 10^4 \rightarrow \log_{10}(10,000) = 4$
- $0.0001 = 10^{-4} \rightarrow -\log_{10}(0.0001) = 4$
- If $[\text{H}^+] = 1 \times 10^{-4} \text{ M}$, then $\text{pH} = 4$


pH & H^+ concentrations in water at 25°C



Why does pH matter?



Minimum pH levels



pH Proficiency Testing



Wastewater manufacturing range: 5.0 to 10.0 units



Wastewater acceptance limits: ± 0.2 units



[NELAC FoPT Tables | PT Ranges and Limits | Waters ERA \(eraqc.com\)](https://eraqc.com)

eraqc.com → Support → NELAC FoPT tables



Example:

DMR-QA 44

pH made-to: 6.49 units


Acceptance limits: 6.29 – 6.69

Outside this range is evaluated as unacceptable

Past DMR-QA results

DMR-QA 44


2,579 Reported values
135 Unacceptable values



Acceptable Unacceptable

DMR-QA 45

2,631 Reported values
124 Unacceptable values




Acceptable Unacceptable

Calibration buffers

pH 4 buffer

- Cat # 127 (500mL)
- Cat # 128 (6 x 500mL)
- Cat # 129 (1L)



pH 7 buffer

- Cat # 131 (500mL)
- Cat # 132 (6 x 500mL)
- Cat # 133 (1L)



pH 10 buffer

- Cat #135 (500mL)
- Cat #136 (6 x 500mL)
- Cat #137 (1L)




Waters ERA products

Proficiency Testing products

Wastewater Cat # 577

Drinking water Cat # 552



Waters™ | ERA

Quality Control Products

Wastewater Cat # 977

Drinking water Cat # 779




pH analysis

- Analysis
 - Pre-measurement checks
 - Proper measurement techniques
- Calibration
 - Use two calibration points
 - Account for temperature effects
 - Verify results with a third-party standard or QC check
- Proper probe care
 - Correct storage practices
 - Routine cleaning
 - Knowing when to replace the probe

Waters™ | ERA




Before you begin pH analysis



- Check probe condition
 - Confirm the probe was stored in an appropriate solution
 - Inspect for crystal buildup
 - Verify internal solution level is sufficient
- Ensure all solutions are at the same temperature
- Check that the fill hole is open
- Confirm Automatic Temperature Compensation (ATC) is enabled

$$E = E^0 - \frac{2.303RT}{F} \text{pH}$$

Probe close-up




pH measurement best practices

- Calibrate on the day of PT analysis
- Submerge probe past the reference junction
- Stir all samples consistently
- Allow readings to stabilize before recording
- Do not wipe the probe dry
- Gently blot or condition the probe between samples

Waters™ | ERA



Calibration of the pH probe



- Use fresh calibration solution
- Perform a two-point calibration only
 - Always include pH 7
 - Use pH 4 or pH 10 as the second buffer
 - Estimate sample pH with pH paper before calibration
 - Enter temperature-corrected buffer values from the bottle
- Enable Automatic Temperature Compensation (ATC), if available

Calibration checks

- Verify calibration slope is between 92% - 102%
 - Target 96% - 101% for PT testing
- Confirm calibration using QC
 - Use a purchased QC sample
 - Sodium tetraborate (pH 9.2) may be used as a check standard
 - Verify with a second source pH buffer
 - Use an additional buffer within the calibration range



Probe storage




Waters™ | ERA

- Follow the probe manufacturer's instructions
- Do **not** store the probe in DI water
- Do **not** allow the probe to dry out
- Keep the bulb fully immersed in manufacturer-recommended storage solution
- If the probe dries out, rehydrate in storage solution for at least 30 minutes

Probe cleaning

- Between samples
 - Do **not** wipe the probe with a tissue
 - Gently blotting is acceptable
- If crystals are present
 - Empty the probe chamber
 - Do **not** add DI water
 - Rinse with probe solution to condition
 - Refill internal solution to the fill line
- Probe-specific notes
 - Gel-filled probes do **not** have replaceable solution
- Cleaning contamination
 - If slime or buildup is present, soak probe in 0.1M acid for 30 minutes
 - Replace internal solution after soaking (when applicable)



pH probe troubleshooting flow

Start → Check mV at pH 7

Is reading within ± 30 mV of 0?

Yes → Proceed to next check

No → Clean and/or rehydrate probe → Recheck



→ Check mV at pH 4 or pH 10

Has mV changed 160-180 mV from pH 7?

Yes → Probe is performing correctly

No → Clean and/or rehydrate probe → Recheck



→ After cleaning/rehydration

Did performance improve?

Yes → Return probe to service

No → Replace probe



Conclusions



pH is measuring a H^+ in solution on a log scale



Proficiency testing for pH has tight limits



pH measurements

Check probe prior to starting
Temperature matters
Calibrate



Probe care

Store properly
Cleaning may help
Sometimes new is the way to go

A background graphic consisting of a network of interconnected nodes and lines, rendered in shades of blue and white, creating a sense of connectivity and data flow. The nodes are represented by small circles of varying sizes, and the lines are thin, light blue lines connecting these nodes.

Waters™ |  ERA™

THANK YOU!