


Preventing Background Contamination in Low-Level TOC Testing

Background TOC doesn't come from one source. It accumulates across six stages between your water system and your reported result. This guide maps the full pathway so you can isolate the real contributor instead of chasing the wrong variable.

The contamination pathway: water system → reported result

Each stage can add background carbon. Elevated blanks usually trace to one or two of these—not all six.



Highlighted stages (2-4) are the most common and controllable contributors in low-level work.

Environment

Airborne VOCs and CO₂ absorb into exposed water surfaces. Open-bench work in high-traffic labs adds measurable background, especially at sub-100 ppb targets.

Key factor: Hood/covered staging vs. open bench.

Container

Residual organics from manufacturing, packaging materials, or storage vary widely between container types. Not all "low-TOC" claims mean the same thing.

Key factor: Certified ultra-low-TOC vs. general lab-grade.

Closure & Septum

Silicone septa and certain cap liners leach organics on contact—especially with extended soak time. Incompatible closures can add 20-50+ ppb in worst cases.

Key factor: Material compatibility and tightening practice.

Operator Handling

Every second a vial sits uncapped, it collects airborne contamination. Inconsistent technique between operators is a top source of blank variability.

Key factor: Uncapped time and fill-then-cap discipline.

Storage & Hold Time

Longer holds give low-level leachables time to accumulate. Temperature and light can accelerate the effect, especially with plastic containers.

Key factor: Time from fill to analysis; storage conditions.

Instrument Pathway

Autosampler needles, transfer tubing, and previous-sample residue can carry organics into your blank. This looks like a container problem but isn't.

Key factor: System blanks vs. vial blanks to separate sources.


Isolating the source matters. Changing containers won't help if the problem is in the autosampler tubing. Running extra rinses won't help if the septum is leaching. Start by comparing vial blanks to system blanks—that one test separates container/handling effects from instrument effects.

What to Manage First

Not every source carries equal weight. This prioritization reflects where labs doing low-level TOC work typically find the most improvement with the least disruption. Standardize one variable at a time and confirm with blanks.

Relative impact on background TOC

Directional guide for low-level work, not survey data.



Tip: If blanks improve after changing containers or closures, the prior supply was a contributor. If they don't, look downstream at instrument or environment.

Bench-ready checklist

Print this. Post it at the bench. Consistent habits reduce blank variability.

Before you open anything

- Confirm container + closure match the workflow (including autosampler fit).
- Stage under a hood or covered area when feasible.

Fill & cap

- Open only the container you're filling next.
- Seal = fill + cap immediately. Minimize uncapped time.
- Avoid contact with the container lip and inside of the closure.

Hold & store

- Return capped containers to covered staging immediately.
- Don't add "just in case" rinse steps your SOP doesn't require.
- Document container format, closure type, and hold conditions.

Quick troubleshooting: where to look first

SYMPTOM	MOST LIKELY SOURCE	FIRST STEP
Blanks consistently high across all operators	Container, closure, or instrument	Compare vial blanks vs. system blanks to isolate
Blanks vary between operators	Handling technique (uncapped time, fill practice)	Standardize fill/cap SOP; observe and time each step
Blanks fine at start, drift up through the run	Hold time / storage, or autosampler carryover	Shorten hold time; run system blanks between samples
Blanks spike after changing vial lot or supplier	Container cleanliness or packaging residue	Run vial blanks on old vs. new lot side by side

Find your starting point

Use the **Vial Cleanliness Risk Checker** to match your target limit, instrument, and workflow to a container recommendation and micro-SOP fixes. Or share your setup with an analyst and we'll help you isolate the variable.

<https://info.eraqc.com/vial-cleanliness-risk-checker-toc>

